



Dimension in mm [Inch]

Type	Connections	X		~Y	
		mm	Inch	mm	Inch
W	1/32"	Ø 25.4	Ø1"	32	1.25"
(Shown) W	1/16"	Ø 31.75	Ø1.25"	44	1.73"
UW	1/16"	Ø 38.1	Ø1.5"	49	1.93"
UW	1/8"	Ø 44.45	Ø1.75"	57.4	2.26"

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen:
MODIFICATIONS:
PART NO.:

Material
MATERIAL

Dimensionen
DIMENSIONS

Projektion
PROJECTION

Allgemein Toleranzen
GENERAL TOLERANCES

SN EN 22768-1/ ISO 2768

**Air actuated, PFAF version
W & UW- Type valves**
with clouse-mount hardware. General dimensions

Masstab
SCALE
1:1

Gezeichnet
DRAWN
07.04.2009
F.P.

Geprüft / Freigabe
CHECKED/ RELEASE
07.04.2009
FCo

CAD Nr.
CAD NO.
506348.idw

Unterverz.:
SUB DIREC.:
KUNDE

Zeichnungs Nr. / DRAWING NO.

VICI[®] VICI AG International, CH-6214 Schenkon

3 - GE - 5063.48