


SHOWN CZUPK



Part#	Description	Bore		L		X1		X2
		mm	Inch	mm	Inch	mm	Inch	
CUC*	Union 1/4-28 to 1/4-28, 1/16" OD tubing	Ø0.25	Ø.010"	-49.5	-1.95"	Ø9.53	Ø 3/8"	Ø9.78 [Ø.385"]
CU*	Union 1/4-28 to 1/4-28, 1/16" OD tubing	Ø0.50	Ø.020"	-49.5	-1.95"	Ø9.53	Ø 3/8"	Ø9.78 [Ø.385"]
CUM*	Union 1/4-28 to 1/4-28, 1/16" OD tubing	Ø0.75	Ø.030"	-49.5	-1.95"	Ø9.53	Ø 3/8"	Ø9.78 [Ø.385"]
CUL*	Union 1/4-28 to 1/4-28, 1/8" OD tubing	Ø1.50	Ø.060"	-50.8	-2.0"	Ø9.53	Ø 3/8"	Ø9.78 [Ø.385"]
CZUC*	Union 1/4-28 to 10-32, 1/16" OD tubing	Ø0.25	Ø.010"	-54.5	-2.15"	Ø9.53	Ø 3/8"	Ø9.65 [Ø.38"]
CZU*	Union 1/4-28 to 10-32, 1/16" OD tubing	Ø0.50	Ø.020"	-54.5	-2.15"	Ø9.53	Ø 3/8"	Ø9.65 [Ø.38"]
CZUM*	Union 1/4-28 to 10-32, 1/16" OD tubing	Ø0.75	Ø.030"	-54.5	-2.15"	Ø9.53	Ø 3/8"	Ø9.65 [Ø.38"]
CZUL*	Union 1/4-28 to 10-32, 1/8" OD tubing	Ø1.0	Ø.040"	-55.0	-2.17"	Ø9.53	Ø 3/8"	Ø9.65 [Ø.38"]
CU4L*	Union 1/4-28 to 1/2-20, 1/8" to 1/4" OD tubing	Ø1.5	Ø.060"	-55.8	-2.20"	Ø19.05	Ø 3/4"	Ø15.75 [Ø.62"]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material	Dimensional
Production	General Tolerances

Part No. _____

**Low pressure unions
polymeric, 1/4-28 to ***

General dimensions

Scale	Drawing No.	20.02.2008	F.P.
5:1	Revision	20.02.2008	FCo
	CAD No.	506322.idw	

Customer: KUNDE

Ziehungs Nr. / Drawing No. 1 - GE - 5063.22

VICI VICI AG International, CH-6214 Schenkon