



SHOWN ZBRU2.5LFPK

Part#	Description	Bore		L		L2		L3		L4		D		X1		X2		X3	
		mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch
ZBRU1.5FPK	Union 1/16 to 1/32", fingertight	Ø0.25	Ø.010"	-45.7	-1.80"	19.56	.770"	11.94	.470"	4.77	.188"	Ø8.0	Ø.3125"	Ø11.1	Ø.437"	Ø9.65	Ø.380"	Ø6.35	Ø.250"
ZBRU1.5LFPK	Union 1/16 to 1/32", fingertight	Ø0.50	Ø.020"	-45.7	-1.80"	19.56	.770"	11.94	.470"	4.77	.188"	Ø8.0	Ø.3125"	Ø11.1	Ø.437"	Ø9.65	Ø.380"	Ø6.35	Ø.250"
ZBRU1.5TFPK	Union 1/16 to 1/32", fingertight	Ø0.79	Ø1/32"	-45.0	-1.77"	18.80	.740"	11.18	.440"	4.77	.188"	Ø8.0	Ø.3125"	Ø11.1	Ø.437"	Ø9.65	Ø.380"	Ø6.35	Ø.250"
ZBRU1.5PK	Union 1/16 to 1/32", hex	Ø0.25	Ø.010"	-35.0	-1.38"	19.56	.770"	11.94	.470"	4.77	.188"	Ø8.0	Ø.3125"	Ø11.1	Ø.437"	6.35 Hex	1/4" Hex	Ø6.35	Ø.250"
ZBRU1.5LPK	Union 1/16 to 1/32", hex	Ø0.50	Ø.020"	-35.0	-1.38"	19.56	.770"	11.94	.470"	4.77	.188"	Ø8.0	Ø.3125"	Ø11.1	Ø.437"	6.35 Hex	1/4" Hex	Ø6.35	Ø.250"
ZBRU1.5TPK	Union 1/16 to 1/32", hex	Ø0.79	Ø1/32"	-34.3	-1.35"	18.80	.740"	11.18	.440"	4.77	.188"	Ø8.0	Ø.3125"	Ø11.1	Ø.437"	6.35 Hex	1/4" Hex	Ø6.35	Ø.250"
ZBRU2.5FPK	Union 1/8 to 1/32", fingertight	Ø0.25	Ø.010"	-42.6	-1.68"	24.26	.955"	13.46	.530"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	Ø6.35	Ø.250"
ZBRU2.5LFPK	Union 1/8 to 1/32", fingertight	Ø0.50	Ø.020"	-42.6	-1.68"	24.26	.955"	13.46	.530"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	Ø6.35	Ø.250"
ZBRU2.5TFPK	Union 1/8 to 1/32", fingertight	Ø0.79	Ø1/32"	-41.4	-1.63"	22.99	.905"	12.19	.480"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	Ø6.35	Ø.250"
ZBRU21CFPK	Union 1/8 to 1/16", fingertight	Ø0.25	Ø.010"	-56.0	-2.20"	26.29	1.035"	15.49	.610"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	Ø9.65	Ø.380"
ZBRU21FPK	Union 1/8 to 1/16", fingertight	Ø0.75	Ø.030"	-56.0	-2.20"	26.29	1.035"	15.49	.610"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	Ø9.65	Ø.380"
ZBRU21LFPK	Union 1/8 to 1/16", fingertight	Ø1.00	Ø.040"	-56.0	-2.20"	26.29	1.035"	15.49	.610"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	Ø9.65	Ø.380"
ZBRU21TFPK	Union 1/8 to 1/16", fingertight	Ø1.59	Ø1/16"	-54.6	-2.15"	25.02	.985"	14.22	.560"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	Ø9.65	Ø.380"
ZBRU21CPK	Union 1/8 to 1/16", hex	Ø0.25	Ø.010"	-45.0	-1.77"	26.29	1.035"	15.49	.610"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	6.35 Hex	1/4" Hex
ZBRU21MPK	Union 1/8 to 1/16", hex	Ø0.50	Ø.020"	-45.0	-1.77"	26.29	1.035"	15.49	.610"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	6.35 Hex	1/4" Hex
ZBRU21PK	Union 1/8 to 1/16", hex	Ø0.75	Ø.030"	-45.0	-1.77"	26.29	1.035"	15.49	.610"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	6.35 Hex	1/4" Hex
ZBRU21LPK	Union 1/8 to 1/16", hex	Ø1.00	Ø.040"	-45.0	-1.77"	26.29	1.035"	15.49	.610"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	6.35 Hex	1/4" Hex
ZBRU21TPK	Union 1/8 to 1/16", hex	Ø1.59	Ø1/16"	-43.7	-1.72"	25.02	.985"	14.22	.560"	3.18	.125"	Ø9.6	Ø.3750"	11.1	.437"	9.53 Hex	3/8" Hex	6.35 Hex	1/4" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen: a) 21.01.09 FP #568, added information
 Modifikationen: b) 23.06.09 FP #582, changed char. ECN 1 0988

Material: DIMENSIONAL TOLERANCES

Prohibition of production: Alignment Tolerances

General TOLERANCES

Part No.: EN 602268-1 ISO 2768

Bulkhead internal reducing unions, high pressure PEEK

Material: 04.03.2008 F.P.
 Scale: 5:1
 Drawing No.: 23.06.2009 F.C.o.
 CAD No.: 506328.idw
 User: KUNDE
 Zeichnungs Nr. / DRAWING NO.: 1 - GE - 5063.28 b

VICI VICI AG International, CH-6214 Schenkon