



CHEMINERT VALVE	X		Y	
	mm	Inch	mm	Inch
C5-* 4 or 6 port	Ø31.7	Ø1.25"	33	1.3"
C5H-* 8 or 10 port	Ø34.9	Ø1.375"	35.1	1.38"
C75-* 4 or 6 port	Ø31.7	Ø1.25"	38.1	1.5"
C75-* 8 or 10 port	Ø34.9	Ø1.375"	40.2	1.58"

MANUAL CONTROLLER PART NO. I-22679

**DIMENSIONS:**  
58 x 92 x 27 mm  
2.27 x 3.62 x 1.06 inch



<small>THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLICIT OR SOLE, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VALCO. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VALCO AND SHALL BE RETURNED UPON DEMAND.</small>		<small>Model: MTEVAL</small> <small>Production: ORIGINAL, TOLERANCES</small>
<small>Author: a) 18.06.88 FP #96, chg. bracket part No.</small> <small>Modifications: b) 12.11.09 FP #652, added part number</small> <small>c) 25.03.13 FP #716, chg. Man. controller design</small>	<small>Project: 19.11.1997 F.P.</small> <small>Design: 25.03.2013 FCo</small> <small>Drawn: 506012.DWG</small>	<small>Scale: 1:1</small> <small>Material: KUNDE</small> <small>Order No.: 1-GE-5060.12 c</small>
<b>Micro electric actuator Multiposition valve</b> <small>with standoff. General dimensions</small>		