





**C62-6186I shown, 6 port 2 pos. valve**



**Dimension in mm [Inch]**

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.				Material MATERIAL	Dimension DIMENSIONS
Änderungen: MODIFICATIONS:				Projection PROJECTION	Allgemein-Toleranzen GENERAL TOLERANCES
PART NO.:					SN EN 22768-1/ISO 2768
<b>4 or 6 Port 2-pos. valve with 1/4-28 fitting, low pressure</b>				Maßstab SCALE 2:1	Gezeichnet DRAWN 22.07.2010 F.P.
Model C62 and C65, General dimensions				Geprüft / Freigegeben CHECKED / RELEASED 23.07.2010 FCo	
Unterverz.: SUB DIR.:				KUNDE	
VICI AG International, CH-6214 Schenkon				Zeichnungs-Nr. / DRAWING NO. <b>2 - GE - 5063.75</b>	