



C72MU-4676 valve shown

Dimension in mm [Inch]

	Part#	Description
20000 psi liq. Coated stainless stator for use with stainless tubing	C72MU-4674	4 port valve, manual
	C72MU-4674D	4 port replacement valve
	C72MU-4676	6 port valve, manual
	C72MU-4676D	6 port replacement valve
For use with fused silica tubing	C72MFSU-4674	4 port valve, manual
	C72MFSU-4674D	4 port replacement valve
	C72MFSU-4676	6 port valve, manual
	C72MFSU-4676D	6 port replacement valve
For 10000 and 15000 psi liq. version dimension identical.		

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: a) 01.09.15 FP #892, New drawing  
MODIFICATIONS:

PART NO.:  
C72MU-467\*

Material MATERIAL  
Dimensionen DIMENSIONS

Projektion PROJECTION  
Allgemein Toleranzen GENERAL TOLERANCES  
SN EN 22768-1/ ISO 2768

Masstab SCALE 1.5:1  
Gezeichnet DRAWN 01.09.2015 F.P.  
Geprüft / Freigabe CHECKED/ RELEASE 01.09.2015 FCo  
CAD Nr. 506426.idw  
CAD NO.

Unterverz.: KUNDE  
SUB DIREC.:  
Zeichnungs Nr. / DRAWING NO.

**20000 psi UHPLC valves**  
**360 µm fittings, 100 µm bore**  
General dimensions

**VICI** VICI AG International, CH-6214 Schenkon

**3 - GE - 5064.26 a**