


SHOWN CZBUPK



\* Different materials, PK = PEEK, KF = CTFE

Part#	Description	Bore		D		L		L2		L3		X1	X2	X3
		mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch			
CBUC*	Bulkhead union 1/4-28 to 1/4-28, 1/16" OD tubing	Ø0.25	Ø.010"	Ø9.5	Ø0.375"	-49.5	-1.95"	17.53	.69"	9.9	.39"	Ø 1/2"	1/2" Hex x 1/8 TH	Ø9.78 [Ø.385"]
CBU*	Bulkhead union 1/4-28 to 1/4-28, 1/16" OD tubing	Ø0.50	Ø.020"	Ø9.5	Ø0.375"	-49.5	-1.95"	17.53	.69"	9.9	.39"	Ø1/2"	1/2" Hex x 1/8 TH	Ø9.78 [Ø.385"]
CBUM*	Bulkhead union 1/4-28 to 1/4-28, 1/16" OD tubing	Ø0.75	Ø.030"	Ø9.5	Ø0.375"	-49.5	-1.95"	17.53	.69"	9.9	.39"	Ø1/2"	1/2" Hex x 1/8 TH	Ø9.78 [Ø.385"]
CBUL*	Bulkhead union 1/4-28 to 1/4-28, 1/8" OD tubing	Ø1.50	Ø.060"	Ø9.5	Ø0.375"	-50.8	-2.0"	17.53	.69"	9.9	.39"	Ø1/2"	1/2" Hex x 1/8 TH	Ø9.78 [Ø.385"]
CZBUC*	Bulkhead union 1/4-28 to 10-32, 1/16" OD tubing	Ø0.25	Ø.010"	Ø8.0	Ø0.313"	-54.5	-2.15"	19.81	.78"	12.2	.48"	Ø1/2"	1/2" Hex x 3/16 TH	Ø9.65 [Ø.38"]
CZBU*	Bulkhead union 1/4-28 to 10-32, 1/16" OD tubing	Ø0.50	Ø.020"	Ø8.0	Ø0.313"	-54.5	-2.15"	19.81	.78"	12.2	.48"	Ø1/2"	1/2" Hex x 3/16 TH	Ø9.65 [Ø.38"]
CZBUM*	Bulkhead union 1/4-28 to 10-32, 1/16" OD tubing	Ø0.75	Ø.030"	Ø8.0	Ø0.313"	-54.5	-2.15"	19.81	.78"	12.2	.48"	Ø1/2"	1/2" Hex x 3/16 TH	Ø9.65 [Ø.38"]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen / MODIFICATIONS: a) 21.01.09 FP #568, added information

Material:	Material:	Dimension:	Dimension:
Material:	Material:	Dimension:	Dimension:

Part No.:

Scale: 5:1

Drawing No.: 19.02.2008

Revision: 21.01.2009

Customer: KUNDE

Drawing No.: 506321 idw

General dimensions

Zusatzinfo Nr. / DRAWING NO.:

VICI VICI AG International, CH-6214 Schenkon

1 - GE - 5063.21 a