



Assembled Part-Nr. without liner **For ASSY-PN with liners select appropriate size, add designation**

C-NEU.5XFPK	.15
C-NEU.5FPK	.20
	.25
	.30
	.35
	.40

Cheminert Nano liner	mm [Inch] ID	mm [Inch] " X " dim.
C-NL.15	0.15 [.006"]	10.16 [.400"] S
C-NL.20	0.20 [.008"]	27.3 [1.075"] L
C-NL.25	0.25 [.010"]	
C-NL.30	0.30 [.012"]	
C-NL.35	0.35 [.014"]	
C-NL.40	0.40 [.016"]	



Dimension in mm [Inch]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: MODIFICATIONS: Projektion PROJECTION Allgemein Toleranzen GENERAL TOLERANCES
 PART NO.: C-NEU.5*FPK.* Gezeichnet DRAWN 23.04.2012 F.P. Geprüft / Freigabe CHECKED/ RELEASE 23.04.2012 FCo
 SN EN 22768-1/ ISO 2768

Union, ext.: Fingertight

Masstab SCALE 3:1 Gezeichnet DRAWN 23.04.2012 F.P. Geprüft / Freigabe CHECKED/ RELEASE 23.04.2012 FCo
 CAD Nr. CAD NO. 506390.idw

General dimensions Unterverz.: SUB DIREC.: KUNDE

VICI® VICI AG International, CH-6214 Schenkon Zeichnungs Nr. / DRAWING NO. 4 - GE - 5063.90