


SHOWN SFVL



Dimension in mm [Inch]

Fitting size	Bore	Prod No
1/16"	0.50 mm	SFV
1/16"	0.50 mm	SFVO
1/16"	0.75 mm	SFVL
1/16"	0.75 mm	SFVOL

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material  
MATERIAL


Dimensionen  
DIMENSIONS

Aenderungen:  
MODIFICATIONS:

Projektion  
PROJECTION

Allgemein Toleranzen  
GENERAL TOLERANCES

PART NO.:  
SFV\*



SN EN 22768-1/ ISO 2768

# Prime / purge valve 1/16", manual version

General dimensions

Masstab SCALE  1:1	Gezeichnet DRAWN	25.02.2010	F.P.
	Geprüft / Freigabe CHECKED/ RELEASE	25.02.2010	FCo
	CAD Nr. CAD NO.	506361.idw	

Unterverz.:  
SUB DIREC.: KUNDE



VICI AG International, CH-6214 Schenkon

Zeichnungs Nr. / DRAWING NO.

## 4 - GE - 5063.61