



8/32" MOUNTING HOLES

MOTOR DRIVER  
PART# 1-22645  
4 FT LONG

POWER CABLE 24VDC

MOUNTING BRACKET  
PART# 1-22596

DIGITAL INPUT  
PART# 1-22537  
6 FT LONG

MANUAL CONTROLLER  
PART NO. 1-22689

DIMENSIONS:  
58 x 92 x 27 mm  
2.28 x 3.62 x 1.06 inch

MANUAL CONTROLLER CABLE  
PART# 1-22537-01  
3 FT LONG

Serie	Description	Fitting	X	L
C22	4 and 6 port	1/4"-28	ø34.93 mm	40.0 mm
2 pos.			ø1.375"	1.575"
C22Z	4, 6, 8 and 10 port	1/16" ZDV	ø34.93 mm	40.0 mm
2 pos.			ø1.375"	1.575"
C22	8 and 10 port	1/4"-28	ø44.2 mm	43.0 mm
2 pos.			ø1.74"	1.69"

2 pos. 4 port C22 valve shown

CONTROL MODULE

DIMENSIONS:

PART# FOR EH: EH2CA  
PART# FOR EP: EP2CA

83 x 96.5 x 26.5 mm  
3.27 x 3.8 x 1.04 inch

4x HOLES FOR 8/32" SCREWS

Dimension in mm [Inch]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Approved: a) 14.08.10 FP #638, chp. control module Modified: b) 25.03.13 FP #716, chp. Man. controller design	Model: MTEVAL Production: 25.02.2002 F.P. Date: 25.03.2013 FCo 506034.DWG	Dimensional: 25.02.2002 F.P. Original: 25.03.2013 FCo 506034.DWG
Part No.:	Material:	Drawing No.:
<b>Micro electric actuator C22 Serie valve</b> with close-mount hdw. General dimensions	Material: none Drawing: 506034.DWG Date: 25.03.2013 FCo	Drawing No.: KUNDE Drawing No.: 1-GE-5060.34 b

VICI® VICI AG International, CH-6214 Scherikon