



6 Mpos. 7 port C5 valve shown

Serie	Description	Fitting	X	L
C5 Mpos.	4 and 6 port	1/16" ZDV	Ø31.75 mm Ø1.25"	47.9 mm 1.885"
C5 Mpos.	8 and 10 port	1/16" ZDV	Ø34.93 mm Ø1.375"	50.0 mm 1.969"
C75H* Mpos.	4 and 6 port	1/16" ZDV	Ø31.75 mm Ø1.25"	52.9 mm 2.083"
C75H* Mpos.	8 and 10 port	1/16" ZDV	Ø34.93 mm Ø1.375"	55.1 mm 2.170"
C75NX* Mpos.	6 port	1/32" ZDV	Ø31.75 mm Ø1.25"	49.1 mm 1.933"
C75NX* Mpos.	8 and 10 port	1/32" ZDV	Ø31.75 mm Ø1.25"	50.4 mm 1.983"

Dimension in mm [Inch]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen:  
MODIFICATIONS:

PART NO.:

**Air actuated with  
Mpos. cheminert valve**

General dimensions

**VICI** VICI AG International, CH-6214 Schenkon

Masstab SCALE	Gezeichnet DRAWN	25.01.2011	F.P.
	Geprüft / Freigabe CHECKED/ RELEASE	26.01.2011	FCo
	CAD Nr. CAD NO.	506385.idw	

Unterverz.:  
SUB DIREC.: KUNDE

Zeichnungs Nr. / DRAWING NO.

**3 - GE - 5063.85**

SN EN 22768-1/ ISO 2768