



shown C45R-8146EUTA valve, 6 pos.

Dimension in mm [Inch]

1/2-20 UNF threaded fitting details for use with 1/4" OD tubing.

Part#	Description	Bore		L	
		mm	Inch	mm	Inch
C45R-8146EUTA	Valve-plug, 6 pos.	Ø3.18	Ø.125"	46.8	1.84"
C45R-8148EUTA	Valve-plug, 8 pos.	Ø2.8	Ø .11"	45.5	1.79"
C45R-9144EUTA	Valve-plug, 4 pos.	Ø4.57	Ø .18"	46.8	1.84"
C45R-9146EUTA	Valve-plug, 6 pos.	Ø3.94	Ø .155"	46.8	1.84"

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: a) 28.09.15 FP #908, New drawing
MODIFICATIONS:

PART NO.: C45R-*

C45R Valve, 1/2-20 UNF Universal electric actuator

General dimensions

Masstab SCALE: 1:1

Gezeichnet DRAWN: 28.09.2015 F.P.
Geprüft / Freigabe CHECKED/ RELEASE: 28.09.2015 FCo
CAD Nr. CAD NO.: 506428.idw

Unterverz.: KUNDE
SUB DIREC.:

Zeichnungs Nr. / DRAWING NO.: 3 - GE - 5064.28 a

VICI® VICI AG International, CH-6214 Schenkon