



SHOWN EZRU.51



SHOWN EZBRU.51



| Part# | Description | Bore | | L | | L2 | | L3 | | D | | X1 | X2 | X3 | X4 |
|-----------|--|-------|--------|-------|--------|-------|-------|------|-------|------|---------|-------|----------|-----------|----------|
| | | mm | Inch | mm | Inch | mm | Inch | mm | Inch | mm | Inch | | | | |
| EZRU.51 | Internal / external reducing union 1/16" to 1/32" | Ø0.25 | Ø0.10" | -25.2 | ~.99" | | | | | | | Ø6.35 | 1/4" Hex | | |
| EZRU.51L | Internal / external reducing union 1/16" to 1/32" | Ø0.50 | Ø0.20" | -25.2 | ~.99" | | | | | | | Ø6.35 | 1/4" Hex | | |
| EZRU.51T | Internal / external reducing union 1/16" to 1/32" | Ø0.89 | Ø0.35" | -24.4 | ~.96" | | | | | | | Ø6.35 | 1/4" Hex | | |
| EZBRU.51 | Bulkhead internal / external reducing union 1/16" to 1/32" | Ø0.25 | Ø0.10" | -26.0 | ~1.02" | 10.29 | .405" | 5.72 | .225" | Ø8.0 | Ø.3125" | Ø6.35 | 1/4" Hex | 7/16" Hex | 1/2" Hex |
| EZBRU.51L | Bulkhead internal / external reducing union 1/16" to 1/32" | Ø0.50 | Ø0.20" | -26.0 | ~1.02" | 10.29 | .405" | 5.72 | .225" | Ø8.0 | Ø.3125" | Ø6.35 | 1/4" Hex | 7/16" Hex | 1/2" Hex |
| EZBRU.51T | Bulkhead internal / external reducing union 1/16" to 1/32" | Ø0.89 | Ø0.35" | -28.7 | ~1.13" | 12.45 | .490" | 7.87 | .310" | Ø8.0 | Ø.3125" | Ø6.35 | 1/4" Hex | 7/16" Hex | 1/2" Hex |

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

| | | |
|--|------------------------|--------------------------|
| Änderungen / MODIFICATIONS: a) 20-01-09 FP #568, added information | | Material: Dimensions: |
| Part No.: | Position / PRODUCTION: | Material: Dimensions: |
| Scale: 5:1 | Date: 13.03.2006 | F.P.: |
| CAD No.: | Date: 20.01.2009 | F.C.: |
| Drawing No.: 506076.idw | Customer: KUNDE | Drawing No.: |
| General dimensions | | Drawing No.: |
| VICI® VICI AG International, CH-6214 Schenkon | | 1 - GE - 5060.76 a |