


Shown FP2Z3M22

* different materials, standard is 300 series stainless steel

Part#	Description	Bore 1	Bore 2	D	L1	L2	L3	X
FP1Z3M21*	One 1/8-27 NPT female pipe to three 1/16" ZDV fittings	Ø1.0 [Ø.040"]	Ø2.0 [Ø.080"]	Ø31.75 [Ø1.25"]	27.3 [1.075"]	14.6 [.575"]	~45.2 [~1.78"]	1/4" Hex
FP1Z3M22*	One 1/8-27 NPT female pipe to three 1/8" ZDV fittings	Ø1.0 [Ø.040"]	Ø2.0 [Ø.080"]	Ø38.1 [Ø1.5"]	27.3 [1.075"]	14.6 [.575"]	~57.1 [~2.25"]	3/8" Hex
FP1Z3M22L*	One 1/8-27 NPT female pipe to three 1/8" ZDV fittings	Ø2.0 [Ø.080"]	Ø2.0 [Ø.080"]	Ø38.1 [Ø1.5"]	27.3 [1.075"]	14.6 [.575"]	~57.1 [~2.25"]	3/8" Hex
FP1Z3M24*	One 1/8-27 NPT female pipe to three 1/4" ZDV fittings	Ø4.6 [Ø.180"]	Ø4.6 [Ø.180"]	Ø44.45 [Ø1.75"]	28.7 [1.13"]	17.3 [.68"]	~66.3 [~2.61"]	1/2" Hex
FP1Z3M41*	One 1/4-18 NPT female pipe to three 1/16" ZDV fittings	Ø1.0 [Ø.040"]	Ø2.0 [Ø.080"]	Ø31.75 [Ø1.25"]	38.6 [1.52"]	19.6 [.77"]	~45.2 [~1.78"]	1/4" Hex
FP1Z3M42*	One 1/4-18 NPT female pipe to three 1/8" ZDV fittings	Ø2.0 [Ø.080"]	Ø2.0 [Ø.080"]	Ø38.1 [Ø1.5"]	38.6 [1.52"]	19.6 [.77"]	~57.1 [~2.25"]	3/8" Hex
FP1Z3M44*	One 1/4-18 NPT female pipe to three 1/4" ZDV fittings	Ø4.6 [Ø.180"]	Ø4.6 [Ø.180"]	Ø44.45 [Ø1.75"]	38.6 [1.52"]	27.2 [1.07"]	~67.6 [~2.66"]	1/2" Hex
FP2Z3M21*	Two 1/8-27 NPT female pipe to three 1/16" ZDV fittings	Ø1.0 [Ø.040"]	Ø2.0 [Ø.080"]	Ø31.75 [Ø1.25"]	27.3 [1.075"]	14.6 [.575"]	~45.2 [~1.78"]	1/4" Hex
FP2Z3M22*	Two 1/8-27 NPT female pipe to three 1/8" ZDV fittings	Ø2.0 [Ø.080"]	Ø2.0 [Ø.080"]	Ø38.1 [Ø1.5"]	27.3 [1.075"]	14.6 [.575"]	~57.1 [~2.25"]	3/8" Hex
FP2Z3M24*	Two 1/8-27 NPT female pipe to three 1/4" ZDV fittings	Ø4.6 [Ø.180"]	Ø4.6 [Ø.180"]	Ø44.45 [Ø1.75"]	28.7 [1.13"]	17.3 [.68"]	~67.6 [~2.66"]	1/2" Hex
FP2Z3M41*	Two 1/4-18 NPT female pipe to three 1/16" ZDV fittings	Ø1.0 [Ø.040"]	Ø2.0 [Ø.080"]	Ø31.75 [Ø1.25"]	38.6 [1.52"]	19.6 [.77"]	~45.2 [~1.78"]	1/4" Hex
FP2Z3M42*	Two 1/4-18 NPT female pipe to three 1/8" ZDV fittings	Ø2.0 [Ø.080"]	Ø4.6 [Ø.180"]	Ø38.1 [Ø1.5"]	38.6 [1.52"]	19.6 [.77"]	~57.1 [~2.25"]	3/8" Hex
FP2Z3M44*	Two 1/4-18 NPT female pipe to three 1/4" ZDV fittings	Ø4.6 [Ø.180"]	Ø4.6 [Ø.180"]	Ø44.45 [Ø1.75"]	38.6 [1.52"]	27.2 [1.07"]	~66.3 [~2.61"]	1/2" Hex



THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material:
 Dimensions:
 Projection:
 General Tolerances:
 SN EN 22768-1 / ISO 2768

Manifold pipe adapters

Scale: 3:1

506082.idw

KUNDE

Zichrungs Nr. / DRAWING NO.
 2 - GE - 5060.82 a

VICI® VICI AG International, CH-6214 Schenkon