

SHOWN ZAOR11

Dimension in mm [Inch]

Part#	Description	Bore		L		X1	X2
		mm	Inch	mm	Inch		
ZAOR11	Internal adaptor 1/16" to 10-32 thread	Ø1.0	Ø.040"	~21.5	~.845"	1/4" Hex	3/8" Hex
ZAOR11HC	Internal adaptor 1/16" to 10-32 thread	Ø1.0	Ø.040"	~21.5	~.845"	1/4" Hex	3/8", 2 Flats OD: Ø0.5"
ZAOR11T	Internal adaptor 1/16" to 10-32 thread	Ø1.65	Ø.065"	~21.5	~.845"	1/4" Hex	3/8" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material
MATERIAL


Dimensionen
DIMENSIONS

Aenderungen: a) 30.08.10 FP #636, added HC material
MODIFICATIONS:

Projektion
PROJECTION

Allgemein Toleranzen
GENERAL TOLERANCES

PART NO.:



SN EN 22768-1/ ISO 2768

Internal adaptor 1/16" to 10-32 thread

General dimensions

Masstab
SCALE

Gezeichnet
DRAWN

22.06.2006

F.P.

3:1

Geprüft / Freigabe
CHECKED/ RELEASE

30.08.2010

FCo

CAD Nr.
CAD NO.

506085.idw

Unterverz.:
SUB DIREC.:

KUNDE

Zeichnungs Nr. / DRAWING NO.

4 - GE - 5060.85 a



VICI AG International, CH-6214 Schenkon