



Dimension in mm [Inch]

Type	Connections	X		~Y		THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.	Material MATERIAL	Dimensionen DIMENSIONS
		mm	Inch	mm	Inch			
W	1/32"	ϕ 25.4	ϕ 1"	32	1.25"	Aenderungen: MODIFICATIONS: a) 23.08.05 FP #461, Redrawn in Inventor	Projektion PROJECTION	Allgemein Toleranzen GENERAL TOLERANCES SN EN 22768-1/ ISO 2768
(Shown) W	1/16"	ϕ 31.75	ϕ 1.25"	44	1.73"			
UW	1/16"	ϕ 38.1	ϕ 1.5"	49	1.93"	PART NO.:	Masstab SCALE 1:1	Gezeichnet DRAWN 03.11.1997 F.P. Geprüft / Freigabe CHECKED/ RELEASE 23.08.2005 FCo CAD Nr. CAD NO. 506009.idw
UW	1/8"	ϕ 44.45	ϕ 1.75"	57.4	2.26"			
Micro electric actuator 2- pos. with standoff. General dimensions							Unterverz.: SUB DIREC.: KUNDE	
VICI AG International, CH-6214 Schenkon							Zeichnungs Nr. / DRAWING NO. 3 - GE - 5060.09 a	