





	Valve Type	Connections	No. of Pos.	X		Y	
				mm	Inch	mm	Inch
 <p>ONE OFFSET PORT</p> <p>ONE ROW OF PORTS</p>	CSD*UW	1/16"	4, 6, 8, 10	ø38.1	ø1.5	~34	~1.34
	CSD*UW	1/16"	12, 16	ø34.9	ø1.375	~32.0	~1.25
	SD*UW	1/8"	4, 6, 8	ø44.45	ø1.75	~38.0	~1.50
 <p>TWO OFFSET PORTS</p> <p>ONE ROW OF PORTS</p>	CST*UW	1/16"	4, 6	ø34.9	ø1.375	~31.5	~1.24
	ST4UW	1/8"	4	ø44.45	ø1.75	~39.5	~1.56



Dimension in mm [Inch]

<small>THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VEGA. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VEGA AND SHALL BE RETURNED UPON DEMAND.</small>		<small>Material:</small> <small>PROJECTION:</small>	<small>Dimension:</small> <small>DIMENSIONS:</small>
<small>Änderungen:</small> a) 07.09.09 FP #592, chg. table data Y		<small>Project:</small>	<small>Approval:</small>
<small>PART NO.:</small>	<small>Scale:</small> 1:5:1	<small>Drawn:</small> 14.06.1999 F.P.	<small>Checked:</small> 07.09.2009 FCo
Micro electric actuator Multiposition Valve <small>with close-mount hdw. General dimensions</small>		<small>Customer:</small> KUNDE	<small>Drawing No. / DRAWING NO.:</small> 1-GE-5060.26 a
 <small>VICI AG International, CH-6214 Schenkon</small>		<small>508026.DWG</small>	