


SCHEMATIC FLOW DIAGRAM

Technical Note 823

CHEMINERT® MODEL C35ZT TRAPPING STREAM SELECTORS

Cheminert Model C35ZT valves are used for multi-column, multi-sample, or multi-trap applications. A typical application, shown below, is the collection of fractions for analysis at a later time.

FIGURE 1: Valve viewed from end opposite handle or actuator



CAUTION: The number of positions indexed by the actuator is half the quantity of numbered ports on the cap:

20 numbered ports + 1 outlet + 1 inlet = 10 positions

24 numbered ports + 1 outlet + 1 inlet = 12 positions

26 numbered ports + 1 outlet + 1 inlet = 13 positions